

Random Hexamer

Random Hexamer; Catalog number: S-6013; **Size:** 50 μl 100 pmoles/ μl ; lyophilized

Notes: Research Use Only.

This product insert declares that this product has been analysed and passed the quality control tests at the time of manufacture.

Shipping: ambient temperature

Store at -20°C after preparation

Random Hexamer Primer is a mixture of single-stranded random hexanucleotides. The sequence is 5'-NNNNNN-3', where N is dA, dT, dG, or dC. The primer is supplied as lyophilized, and is enough for 50 first strand cDNA synthesis reactions.

Preparation before first use:

Centrifuge the tube at 13000 rpm for 2 minutes. Add the volume indicated on the tube from nuclease-free water, vortex for 1 min, and let it seat for 3 minutes. Mix by vortexing for 1 min and briefly spin. This gives a 100 μM (100 pmoles/ μl) solution.

Guidelines for reverse transcription reaction

1. Thaw RNA and other materials on ice.
2. Add the reagents in the following order:

| | |
|---|--------------------------|
| • RNA (100 nanogram to 5 microgram) | ... μl |
| • Random hexamer primer; 1 μl (200 ng) | 1 μl (200 ng) |
| • H ₂ O | up to 13 μl |
3. Incubate the tube at 70°C for 5 min.
4. Quickly chill the tube on ice.
5. Add the followings in order:

| | |
|--------------------------|-----------------|
| • RT buffer (5x) | 4 μl |
| • dNTPs mix (10 mM each) | 2 μl |
6. Vortex, spin and incubate at 37°C for 5 min.
7. Add 1 μl (200 U) of RT enzyme.
8. Vortex, spin and incubate at the indicated temperature (for MMLV, 42°C) for 60 min. cDNA is made at the end of this period.
9. Stop the reaction by heating at 70°C for 10 min.
10. Chill on ice.
11. Store the cDNA at -20°C, or quickly set up a PCR reaction using cDNA as the template.