

DENazist dNTP mix

Shipping: blue ice
Store at -20°C



Cat. No.: S-6010



معرفی و کاربردها: مخلوط بازهای DNA (شامل بر dGTP، dATP، dTTP، dCTP و dGTP) در واکنش های مختلف بیولوژی مولکولی از جمله PCR و RT مورد استفاده قرار می گیرد. غلظت این بازها در این مخلوط هر کدام ۱۰ میلی مولار می باشد. بعضی موارد استفاده عبارتند از:

- PCR and qPCR
- Reverse transcription
- DNA labeling

Guide to setup a 25 ul PCR reaction

- | | |
|-------------------------|------------------------|
| ■ DNA template | as required |
| ■ DNA polymerase | 0.5 to 5 U |
| ■ DNA polymerase buffer | final 1x concentration |
| ■ Primers (10 uM each) | 1 ul |
| ■ dNTP mix (10 mM each) | 0.5 ul |

DENazist dNTP mix

Cat. No.: S-6010



Shipping: blue ice
Store at -20°C

Description:

DENazist dNTP mix (including dATP, dTTP, dCTP, and dGTP) is used many molecular biology applications including PCR and RT. The concentration of each base is 10 mM. Some applications include:

- PCR and qPCR
- Reverse transcription
- DNA labeling

Guide to setup a 25 ul PCR reaction

- | | |
|-------------------------|------------------------|
| ■ DNA template | as required |
| ■ DNA polymerase | 0.5 to 5 U |
| ■ DNA polymerase buffer | final 1x concentration |
| ■ Primers (10 uM each) | 1 ul |
| ■ dNTP mix (10 mM each) | 0.5 ul |

Notes: Research Use Only

This product insert declares that this products has been analyzed and passed the quality control tests at the time of manufacture.